

Product datasheet

VPS45A peptide (294-307) ab45784

製品の概要

製品名 VPS45A peptide (294-307)

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 C-FQKKKPKEQQKLES

領域 294 to 307

特性

Our [Abpromise guarantee](#) covers the use of **ab45784** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-VPS45A antibody ([ab40853](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

関連情報

機能 May play a role in vesicle-mediated protein trafficking from the Golgi stack through the trans-Golgi network.

組織特異性 Ubiquitous. Expression was highest in testis, heart and brain, intermediate in kidney, spleen, prostate, ovary, small intestine and thymus and low in lung, skeletal muscle, placenta, colon, pancreas, peripheral blood leukocytes and liver.

配列類似性 Belongs to the STXBP/unc-18/SEC1 family.

細胞内局在 Golgi apparatus membrane. Endosome membrane. Associated with Golgi/endosomal vesicles and the trans-Golgi network.

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors