

Product datasheet

Anti-XIRP2 antibody ab122234

画像数 1

製品の概要

製品名	Anti-XIRP2 antibody
製品の詳細	Rabbit polyclonal to XIRP2
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	antigen sequence: HAPPTYEDVI AGHILDISDS PKEVRKNFQK TWQESGRVFK GLGYATADAS ATEMRTTFQE ESAFISEAAA PRQGNMYTLS KDSL SNGVP, corresponding to amino acids 3278-3366 of Human XIRP2 (UniProt ID: A4UGR9). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human kidney tissue

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122234** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

機能	Protects actin filaments from depolymerization.
配列類似性	Belongs to the Xin family. Contains 28 Xin repeats.
ドメイン	Xin repeats bind F-actin.
細胞内局在	Cell junction. Colocalizes with actin stress fibers.

画像



ab122234, at a dilution of 1/75, staining XIRP2 in paraffin embedded Human kidney tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-XIRP2 antibody (ab122234)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors