

Product datasheet

Anti-XIRP2 antibody ab122234

画像数 1

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-XIRP2 antibody |
| 製品の詳細 | Rabbit polyclonal to XIRP2 |
| アプリケーション | 適用あり: IHC-P |
| 種交差性 | 交差種: Human |
| 免疫原 | antigen sequence: HAPPTYEDVI AGHILDISDS PKEVRKNFQK TWQESGRVFK GLGYATADAS ATEMRTTFQE ESAFISEAAA PRQGNMYTLS KDLSLNGVP, corresponding to amino acids 3278-3366 of Human XIRP2 (UniProt ID: A4UGR9). Run BLAST with ExPASy Run BLAST with NCBI |
| ポジティブ・コントロール | Human kidney tissue |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| バッファー | pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122234** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| IHC-P | | 1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. |

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|------|
|----------|-----------|------|

ターゲット情報

| | |
|-------|--|
| 機能 | Protects actin filaments from depolymerization. |
| 配列類似性 | Belongs to the Xin family. Contains 28 Xin repeats. |
| ドメイン | Xin repeats bind F-actin. |
| 細胞内局在 | Cell junction. Colocalizes with actin stress fibers. |

画像



ab122234, at a dilution of 1/75, staining XIRP2 in paraffin embedded Human kidney tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-XIRP2 antibody (ab122234)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors