

## Product datasheet

# Anti-ZNF776 antibody ab67522

### 画像数 1

#### 製品の概要

製品名	Anti-ZNF776 antibody
製品の詳細	Mouse polyclonal to ZNF776
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein (Human)
ポジティブ・コントロール	ZNF776 transfected lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab67522** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 50 kDa (predicted molecular weight: 60 kDa).

#### ターゲット情報

機能	May be involved in transcriptional regulation.
配列類似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 11 C2H2-type zinc fingers. Contains 1 KRAB domain.
細胞内局在	Nucleus.

## 画像



Western blot - ZNF776 antibody (ab67522)

**All lanes :** Anti-ZNF776 antibody (ab67522)  
at 1/500 dilution

**Lane 1 :** ZNF776 transfected lysate

**Lane 2 :** Non-transfected lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat anti-mouse IgG(H&L)-HRP  
conjugate at 1/2500 dilution

**Predicted band size:** 60 kDa

**Observed band size:** 50 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors