

Product datasheet

Anti-Interferon alpha 1 antibody [AE3] ab10077

★★★★☆ 1 Abreviews

製品の概要

製品名	Anti-Interferon alpha antibody [AE3]
製品の詳細	Mouse monoclonal [AE3] to Interferon alpha
由来種	Mouse
アプリケーション	適用あり: WB, Sandwich ELISA, ELISA
種交差性	交差種: Human
免疫原	Recombinant full length protein (Human).
特記事項	<p>Concentration varies from lot to lot and can be provided on request.</p> <p>Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. If you would like help finding an alternative product, please do not hesitate to contact our scientific support team.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.2
精製度	Protein A purified
特記事項(精製)	Purity tested by electrophoresis.
ポリ/モノ	モノクローナル
クローン名	AE3
ミエローマ	x63-Ag8.653
アイソタイプ	IgG2a

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab10077** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB	★★★★☆	Use at an assay dependent dilution. Predicted molecular weight: 23 kDa.
Sandwich ELISA		Use a concentration of 5 µg/ml. Can be paired for Sandwich ELISA with Mouse monoclonal [AB5] to Interferon alpha (HRP) (ab10078) . Can be used as Capture antibody with recommended pair. Detection limit 5pg/ml.
ELISA		Use at an assay dependent dilution.

ターゲット情報

機能	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
配列類似性	Belongs to the alpha/beta interferon family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors