

## Product datasheet

# Anti-TFIIB antibody ab101280

### 画像数 1

#### 製品の概要

製品名	Anti-TFIIB antibody
製品の詳細	Rabbit polyclonal to TFIIB
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IP <b>適用なし:</b> WB
種交差性	<b>交差種:</b> Human <b>交差が予測される動物種:</b> Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Dog, Turkey, Pig, Xenopus laevis, Chimpanzee, Rhesus monkey, Gorilla, Orangutan, Xenopus tropicalis
免疫原	A region within amino acids 100-150 of Human TFIIB (NP_001505.1)
ポジティブ・コントロール	Jurkat whole cell lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium Azide Constituents: Tris citrate/phosphate buffer, pH 7-8
精製度	Immunogen affinity purified
特記事項 (精製)	ab101280 was affinity purified using an epitope specific to TFIIB immobilized on solid support.
ポリモノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab101280** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IP		Use at 2-5 µg/mg of lysate.


## 追加情報

Is unsuitable for WB.

## ターゲット情報

機能	General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.
配列類似性	Belongs to the TFIIIB family. Contains 1 TFIIIB-type zinc finger.
細胞内局在	Nucleus.

## 画像



Detection of TFIIIB by Western Blot of Immunoprecipitate.

anti TFIIIB, at 1  $\mu\text{g/ml}$ , staining TFIIIB in Jurkat whole cell lysate immunoprecipitated using ab101280 at 3  $\mu\text{g/mg}$  lysate (1 mg/IP; 20% of IP loaded/lane). Detection:

Chemiluminescence with an exposure time of 3 seconds.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors