

Product datasheet

Anti-TFIIB antibody ab101280

画像数 1

製品の概要

製品名	Anti-TFIIB antibody
製品の詳細	Rabbit polyclonal to TFIIB
アプリケーション	適用あり: IP 適用なし: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Dog, Turkey, Pig, Xenopus laevis, Chimpanzee, Rhesus monkey, Gorilla, Orangutan, Xenopus tropicalis
免疫原	A region within amino acids 100-150 of Human TFIIB (NP_001505.1)
ポジティブ・コントロール	Jurkat whole cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium Azide Constituents: Tris citrate/phosphate buffer, pH 7-8
精製度	Immunogen affinity purified
特記事項 (精製)	ab101280 was affinity purified using an epitope specific to TFIIB immobilized on solid support.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab101280** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
IP		Use at 2-5 µg/mg of lysate.

追加情報 Is unsuitable for WB.

ターゲット情報

機能	General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.
配列類似性	Belongs to the TFIIIB family. Contains 1 TFIIIB-type zinc finger.
細胞内局在	Nucleus.

画像



Detection of TFIIIB by Western Blot of Immunoprecipitate.
anti TFIIIB, at 1 μ g/ml, staining TFIIIB in Jurkat whole cell lysate immunoprecipitated using ab101280 at 3 μ g/mg lysate (1 mg/IP; 20% of IP loaded/lane). Detection: Chemiluminescence with an exposure time of 3 seconds.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors