

## Product datasheet

# Anti-GM130 antibody ab103416

★★★★☆ 1 Abreviews [画像数 2](#)

### 製品の概要

製品名	Anti-GM130 antibody
製品の詳細	Rabbit polyclonal to GM130
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB, IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 832-862 of Human GM130.
ポジティブ・コントロール	HepG2 cell lysate. Human hepatocarcinoma.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab103416 is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab103416** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB	★★★★☆	1/100 - 1/500. Predicted molecular weight: 112 kDa.

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## ターゲット情報


機能	Golgi auto-antigen; probably involved in maintaining cis-Golgi structure.
配列類似性	Belongs to the GOLGA2 family.
ドメイン	Extended rod-like protein with coiled-coil domains.
細胞内局在	Golgi apparatus > Golgi stack membrane.

## 画像



Anti-GM130 antibody (ab103416) at 1/100 dilution + HepG2 cell lysate at 35 µg

**Predicted band size:** 112 kDa



ab103416, at 1/50 dilution, staining GM130 in formalin-fixed, paraffin-embedded Human hepatocarcinoma by immunohistochemistry using a peroxidase conjugated secondary antibody and DAB staining.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors