

Product datasheet

Anti-ICK antibody ab103817

画像数 2

製品の概要

製品名	Anti-ICK antibody
製品の詳細	Rabbit polyclonal to ICK
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1-632 of Human ICK (NP_055735.1)
ポジティブ・コントロール	Human placenta tissue lysate and ICK transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab103817** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500. Predicted molecular weight: 71 kDa.

ターゲット情報


関連性 ICK may play a key role in the development of multiple organ systems and particularly in cardiac

development.

細胞内局在

Cytoplasmic and Nuclear


画像



Western blot - ICK antibody (ab103817)

Anti-ICK antibody (ab103817) at 1/500 dilution + Human placenta tissue lysate at 50 μ g

Predicted band size: 71 kDa



Western blot - ICK antibody (ab103817)

All lanes : Anti-ICK antibody (ab103817) at 1/500 dilution

Lane 1 : ICK transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 μ g per lane.

Predicted band size: 71 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors