

## Product datasheet

# Anti-DAP Kinase 1 antibody [DAPK-55] ab10443

## 2 References

### 製品の概要

製品名	Anti-DAP Kinase 1 antibody [DAPK-55]
製品の詳細	Mouse monoclonal [DAPK-55] to DAP Kinase 1
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human, Monkey
免疫原	Recombinant full length protein.
ポジティブ・コントロール	COS-7 whole cell extract.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 15mM Sodium Azide Constituents: Ascites
精製度	Ascites
ポリ/モノ	モノクローナル
クローン名	DAPK-55
アイソタイプ	IgG1

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab10443** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/250. Detects a band of approximately 160 kDa.

<b>機能</b>	Calcium/calmodulin-dependent serine/threonine kinase which acts as a positive regulator of apoptosis.
<b>配列類似性</b>	Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. DAP kinase subfamily. Contains 10 ANK repeats. Contains 1 death domain. Contains 1 protein kinase domain.
<b>翻訳後修飾</b>	Ubiquitinated by the BCR(KLHL20) E3 ubiquitin ligase complex, leading to its degradation by the proteasome.
<b>細胞内局在</b>	Cytoplasm. Colocalizes with the actin filament system.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors