

## Product datasheet

# Anti-SCRN2 antibody ab105356

画像数 3

### 製品の概要

製品名	Anti-SCRN2 antibody
製品の詳細	Rabbit polyclonal to SCRN2
アプリケーション	適用あり: ICC/IF, WB, ICC
種交差性	交差種: Mouse, Rat, Human
免疫原	An 18 amino acid synthetic peptide, from a region within the C-terminus of Human SCRN2 (EAW94970).
ポジティブ・コントロール	293 cells

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (stable for up to 12 months). Store at -20°C or -80°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105356** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 47 kDa.
ICC		Use a concentration of 20 µg/ml.

## ターゲット情報

### 配列類似性


Belongs to the peptidase C69 family. Secernin subfamily.

## Anti-SCRN2 antibody 画像



Immunocytochemistry/ Immunofluorescence - Anti-SCRN2 antibody (ab105356)

Immunofluorescence of SCRN2 in 293 cells with ab105356 at 20 µg/mL.



Western blot - SCRN2 antibody (ab105356)

**Lane 1** : Anti-SCRN2 antibody (ab105356) at 0.5 µg/ml

**Lane 2** : Anti-SCRN2 antibody (ab105356) at 1 µg/ml

**Lane 1** : 293 cell lysate

**Lane 2** : 293 cell lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size** : 47 kDa



Immunocytochemistry - SCRN2 antibody (ab105356)

ab105356, at 20 µg/ml, staining SCRN2 in 293 cells by immunocytochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors