

Product datasheet

Anti-SCRN2 antibody ab105356

画像数 3

製品の概要

| | |
|--------------|--|
| 製品名 | Anti-SCRN2 antibody |
| 製品の詳細 | Rabbit polyclonal to SCRN2 |
| アプリケーション | 適用あり: ICC/IF, WB, ICC |
| 種交差性 | 交差種: Mouse, Rat, Human |
| 免疫原 | An 18 amino acid synthetic peptide, from a region within the C-terminus of Human SCRN2 (EAW94970). |
| ポジティブ・コントロール | 293 cells |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at 4°C (stable for up to 12 months). Store at -20°C or -80°C. |
| バッファー | Preservative: 0.02% Sodium Azide Constituents: PBS |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105356** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| ICC/IF | | Use a concentration of 20 µg/ml. |
| WB | | Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 47 kDa. |
| ICC | | Use a concentration of 20 µg/ml. |

ターゲット情報

配列類似性


Belongs to the peptidase C69 family. Secernin subfamily.

画像



Immunofluorescence of SCRN2 in 293 cells with ab105356 at 20 ug/mL.

Immunocytochemistry/ Immunofluorescence - Anti-SCRN2 antibody (ab105356)



Lane 1 : Anti-SCRN2 antibody (ab105356) at 0.5 µg/ml

Lane 2 : Anti-SCRN2 antibody (ab105356) at 1 µg/ml


Lane 1 : 293 cell lysate

Lane 2 : 293 cell lysate

Lysates/proteins at 15 µg per lane.

Western blot - SCRN2 antibody (ab105356)

Predicted band size : 47 kDa



ab105356, at 20ug/ml, staining SCRN2 in 293 cells by immunocytochemistry.

Immunocytochemistry - SCRN2 antibody (ab105356)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors