

## Product datasheet

# Anti-SCRN2 antibody ab105356

### 画像数 3

#### 製品の概要

製品名	Anti-SCRN2 antibody
製品の詳細	Rabbit polyclonal to SCRN2
由来種	Rabbit
アプリケーション	<b>適用あり:</b> ICC/IF, WB, ICC
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	An 18 amino acid synthetic peptide, from a region within the C-terminus of Human SCRN2 (EAW94970).
ポジティブ・コントロール	293 cells

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (stable for up to 12 months). Store at -20°C or -80°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105356** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 47 kDa.
ICC		Use a concentration of 20 µg/ml.

## ターゲット情報

### 配列類似性


Belongs to the peptidase C69 family. Secernin subfamily.

## 画像



Immunofluorescence of SCRIN2 in 293 cells with ab105356 at 20 ug/mL.

Immunocytochemistry/ Immunofluorescence - Anti-SCRIN2 antibody (ab105356)



**Lane 1** : Anti-SCRIN2 antibody (ab105356) at 0.5 µg/ml


**Lane 2** : Anti-SCRIN2 antibody (ab105356) at 1 µg/ml

**All lanes** : 293 cell lysate

Lysates/proteins at 15 µg per lane.

**Predicted band size:** 47 kDa

Western blot - SCRIN2 antibody (ab105356)



ab105356, at 20ug/ml, staining SCRIN2 in 293 cells by immunocytochemistry.

Immunocytochemistry - SCRIN2 antibody (ab105356)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors