

## Product datasheet

# Anti-C18orf25 antibody ab107064

### 画像数 2

#### 製品の概要

製品名	Anti-C18orf25 antibody
製品の詳細	Rabbit polyclonal to C18orf25
アプリケーション	<b>適用あり:</b> WB, IHC-P
種交差性	<b>交差種:</b> Mouse, Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 163-191 of Human C18orf25 (NP_659492.1, NP_001008240.1).
ポジティブ・コントロール	Mouse heart tissue lysate; Human bladder carcinoma tissue.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
特記事項 (精製)	ab107064 was purified through a protein A column, followed by peptide affinity purification.
ポリモノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab107064** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 43 kDa.
IHC-P		1/50 - 1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## ターゲット情報

### 翻訳後修飾


Phosphorylated upon DNA damage, probably by ATM or ATR.

### 画像



Anti-C18orf25 antibody (ab107064) at 1/100 dilution + Mouse heart tissue lysate at 35 µg

**Predicted band size** : 43 kDa



ab107064 at 1/50 dilution staining C18orf25 in Formalin-fixed, Paraffin-embedded Human bladder carcinoma tissue by Immunohistochemistry followed by peroxidase-conjugated secondary antibody and DAB staining.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors