

## Product datasheet

# Anti-NAT14 antibody ab107546

★★★★☆ 2 Abreviews [画像数 1](#)

### 製品の概要

製品名	Anti-NAT14 antibody
製品の詳細	Rabbit polyclonal to NAT14
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Human
免疫原	KLH conjugated synthetic peptide designed between amino acids 12-41 from the N-terminal region of Human NAT14 (NP_065111.1)
ポジティブ・コントロール	Mouse bladder tissue.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab107546 is purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab107546** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 22 kDa.

<b>関連性</b>	Belongs to the camello family and contains 1 N-acetyltransferase domain. Expressed in K562 cells, HeLa cells and brain. Probable acetyltransferase that binds the 5'-GGACTACAG-3' sequence of coproporphyrinogen oxidase promoter. Able to activate transcription of a reporter construct in vitro.
<b>細胞内局在</b>	Membrane; Single pass membrane protein

## 画像



Anti-NAT14 antibody (ab107546) at 1/100 dilution + mouse bladder tissue lysates at 35 µg

**Predicted band size:** 22 kDa

Western blot - NAT14 antibody (ab107546)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors