

Product datasheet

Anti-CD8B antibody [CT-CD8b] (FITC) ab112168

製品の概要

製品名	Anti-CD8B antibody [CT-CD8b] (FITC)
製品の詳細	Rat monoclonal [CT-CD8b] to CD8B (FITC)
由来種	Rat
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	Mouse CD8B
ポジティブ・コントロール	Mouse Strain: BALB/c (Thymus)

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.02% Sodium azide Constituents: 99% PBS, 0.5% BSA
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	CT-CD8b
アイソタイプ	IgG2a

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab112168** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 2µg for 10 ⁶ cells. ab18446 - Rat monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

ターゲット情報

細胞内局在

Isoform 1: Cell membrane; Single-pass type I membrane protein (Probable). Isoform 2: Cell membrane; Single-pass type I membrane protein (Probable). Isoform 3: Secreted (Probable). Isoform 4: Cell membrane; Single-pass type I membrane protein (Probable). Isoform 5: Cell membrane; Single-pass type I membrane protein (Probable). Isoform 6: Secreted (Probable). Isoform 7: Secreted (Probable). Isoform 8: Secreted (Probable).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors