

## Product datasheet

# Goat Anti-Cat IgA H&L ab112793

### 製品の概要

製品名	Goat Anti-Cat IgA H&L
由来種	Goat
ターゲット生物種	Cat
アプリケーション	<b>適用あり:</b> WB, ELISA, IHC-P, ICC/IF
免疫原	Other Immunogen Type corresponding to Cat IgA.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.1% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
特記事項 (精製)	Solid phase adsorbed to ensure class specificity. ab112793 was isolated by affinity chromatography using antigen coupled to agarose beads.
ポリモノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab112793** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/30000. Colorimetric: 1/1,000 - 1/10,000 Chemiluminescent: 1/1,000 - 1/30,000
ELISA		1/100 - 1/30000. Coating: 1/100 - 1/500 Primary: 1/1000 - 1/30,000
IHC-P		1/200 - 1/2000.

アプリケーション

Abreviews

特記事項

ICC/IF

1/200 - 1/2000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors