

Product datasheet

Anti-Myeloid Precursor Cells antibody [OX-82] (FITC) ab114076

1 References 画像数 1

製品の概要

製品名	Anti-Myeloid Precursor Cells antibody [OX-82] (FITC)
製品の詳細	Mouse monoclonal [OX-82] to Myeloid Precursor Cells (FITC)
由来種	Mouse
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Rat
免疫原	Anemic Rat bone marrow.
ポジティブ・コントロール	Bone marrow cells.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.02% Sodium azide Constituents: 99% PBS, 0.5% BSA
精製度	Protein G purified
特記事項 (精製)	ab114076 was purified from ascitic fluid via Protein G Chromatography.
ポリモノ	モノクローナル
クローン名	OX-82
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab114076** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

Use 0.1µg for 10<sup>6</sup> cells.

[ab91356](#) - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## ターゲット情報


### 関連性

ER-MP58 is expressed by all bone marrow derived M-CSF- and GM-CSF-responsive myeloid blood cell precursors. The expression of ER-MP58 remains at a high level throughout the precursor/monocyte stage and is down regulated upon maturation into mature macrophages. The ER-MP58 antigen is used to distinguish between early myeloid committed cells and other haematopoietic progenitors cells in the BM.

### 細胞内局在

Cell Membrane

## 画像



ab114076, at 0.1 µg/10<sup>6</sup> cells, staining Rat Myeloid Lineage in Rat (strain: Wistar) bone marrow cells by Flow Cytometry.

Isotypic Control: FITC Mouse IgG1

Percentage of cells stained above control:

21.5%

Flow Cytometry - Anti-Myeloid Precursor Cells antibody [OX-82] (FITC) (ab114076)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors