

Product datasheet

Anti-Hsp90 alpha antibody [9D2] (HRP) ab115649

画像数 1

製品の概要

製品名	Anti-Hsp90 alpha antibody [9D2] (HRP)
製品の詳細	Rat monoclonal [9D2] to Hsp90 alpha (HRP)
由来種	Rat
標識	HRP
アプリケーション	適用あり: WB
種交差性	交差種: Chicken, Human
免疫原	Purified Hsp90 alpha, isolated from Human therapeutic orchiectomy specimens.
ポジティブ・コントロール	Hsp90 alpha native and recombinant proteins; HeLa Cell Lysate (heat shocked).

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: 99% PBS
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	9D2
アイソタイプ	IgG2a

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab115649** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/2000 - 1/50000. Detects a band of approximately 90 kDa (predicted molecular weight: 85 kDa).

## ターゲット情報

機能	Molecular chaperone. Has ATPase activity.
配列類似性	Belongs to the heat shock protein 90 family.
翻訳後修飾	ISGylated.
細胞内局在	Cytoplasm. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

## 画像



Western blot - Anti-Hsp90 alpha antibody [9D2]  
(HRP) (ab115649)

**All lanes :** Anti-Hsp90 alpha antibody [9D2]  
(HRP) (ab115649) at 1/2000 dilution

**Lane 1 :** Hsp90 alpha native protein

**Lane 2 :** Hsp90 alpha recombinant protein

**Lane 3 :** HeLa Cell Lysate (Heat Shocked)

**Predicted band size:** 85 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors