


Product datasheet

Anti-CLIP170 antibody ab116301

画像数 1

製品の概要

製品名	Anti-CLIP170 antibody
製品の詳細	Rabbit polyclonal to CLIP170
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Horse, Guinea pig, Cow, Cat, Dog 
免疫原	Synthetic peptide corresponding to a region within C-terminal amino acids 1064-1113 (LRKQADKAKA AQTAEAMQI MEQMTKEKTE TLASLEDTKQ TNAKLQNELD) of Human CLIP170 (NP_002947). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human Fetal Liver lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 97% PBS, 2% Sucrose
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab116301** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
----------	-----------	------

WB	Use a concentration of 0.125 µg/ml. Predicted molecular weight: 162 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	---

ターゲット情報

機能	Seems to be a intermediate filament associated protein that links endocytic vesicles to microtubules.
組織特異性	Highly expressed in the Reed-Sternberg cells of Hodgkin disease.
配列類似性	Contains 2 CAP-Gly domains.
細胞内局在	Cytoplasm. Cytoplasm > cytoskeleton. Associated with the cytoskeleton. Localizes to the tips of microtubules.

画像



Anti-CLIP170 antibody (ab116301) at 0.125 μ g/ml + Human Fetal Liver lysate at 10 μ g

Predicted band size : 162 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors