

Product datasheet

Anti-KLHL2 antibody [8D8] ab117934

画像数 1

製品の概要

| | |
|--------------|--|
| 製品名 | Anti-KLHL2 antibody [8D8] |
| 製品の詳細 | Mouse monoclonal [8D8] to KLHL2 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Human |
| 免疫原 | Recombinant fragment, corresponding to amino acids 1-100 and 494-593 of Human KLHL2 produced in HEK293T cells (NP_009177). |
| ポジティブ・コントロール | HEK293T cell lysate transfected with pCMV6-ENTRY KLHL2 cDNA. |
| 特記事項 | Stable for 12 months from date of receipt. |

製品の特性

| | |
|-----------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| バッファー | pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 1% BSA, 50% Glycerol |
| 精製度 | Protein A purified |
| 特記事項 (精製) | ab117934 is purified from Mouse ascites fluids by affinity chromatography. |
| ポリ/モノ | モノクローナル |
| クローン名 | 8D8 |
| アイソタイプ | IgG2b |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab117934** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | 1/500. Predicted molecular weight: 65 kDa. |

ターゲット情報

| | |
|--------------|--|
| 機能 | May play a role in organizing the actin cytoskeleton of the brain cells. |
| 組織特異性 | Predominantly expressed in brain. |
| 配列類似性 | Contains 1 BTB (POZ) domain. |
| | Contains 6 Kelch repeats. |
| 細胞内局在 | Cytoplasm > cytoskeleton. |

画像



All lanes : Anti-KLHL2 antibody [8D8] (ab117934) at 1/500 dilution

Lane 1 : HEK293T cell lysate transfected with pCMV6-ENTRY control cDNA

Lane 2 : HEK293T cell lysate transfected with pCMV6-ENTRY KLHL2 cDNA

Lysates/proteins at 5 µg per lane.

Predicted band size: 65 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC205347)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors