

## Product datasheet

# Anti-C10orf82 antibody [1G4] ab117935

1 References [画像数 2](#)

### 製品の概要

製品名	Anti-C10orf82 antibody [1G4]
製品の詳細	Mouse monoclonal [1G4] to C10orf82
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB, ICC/IF
種交差性	<b>交差種:</b> Human
免疫原	Recombinant full length Human C10orf82 produced in HEK293T cells (NP_653262).
ポジティブ・コントロール	HEK293T cell lysate transfected with pCMV6-ENTRY C10orf82; COS7 cells transiently transfected by pCMV6-ENTRY C10orf82.
特記事項	Stable for 12 months from date of receipt.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 1% BSA, 50% Glycerol
精製度	Protein G purified
特記事項 (精製)	ab117935 is purified from Mouse ascites fluids by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	1G4
アイソタイプ	IgG1

### アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab117935** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/2000. Predicted molecular weight: 26 kDa.
ICC/IF		1/100.


## ターゲット情報

## 画像



Immunocytochemistry/ Immunofluorescence - Anti-C10orf82 antibody [1G4] (ab117935)

ab117935 at 1/100 dilution staining C10orf82 in COS7 cells transiently transfected by pCMV6-ENTRY C10orf82 by Immunofluorescence.



Western blot - Anti-C10orf82 antibody [1G4] (ab117935)

**All lanes :** Anti-C10orf82 antibody [1G4] (ab117935) at 1/2000 dilution

**Lane 1 :** HEK293T cell lysate transfected with pCMV6-ENTRY control cDNA

**Lane 2 :** HEK293T cell lysate transfected with pCMV6-ENTRY C10orf82 cDNA

Lysates/proteins at 5 µg per lane.

**Predicted band size:** 26 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC200834)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors