

## Product datasheet

# Anti-IQCH antibody ab118092

画像数 1

### 製品の概要

製品名	Anti-IQCH antibody
製品の詳細	Mouse monoclonal to IQCH
アプリケーション	<b>適用あり:</b> ELISA, Sandwich ELISA
種交差性	<b>交差種:</b> Recombinant fragment <b>交差が予測される動物種:</b> Human
免疫原	Recombinant fragment, corresponding to internal sequence amino acids 301-401 of Human IQCH (NP_073621). Mol Wt 37.92 kDa inclusive of tag. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Recombinant protein

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Constituent: 99% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab118092** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.

## アプリケーション

Sandwich ELISA

## Abreviews


## 特記事項

Use a concentration of 10 µg/ml.

## ターゲット情報

機能	May play a regulatory role in spermatogenesis.
組織特異性	Expressed in fetal and adult testis, in germ cells but not somatic cells.
配列類似性	Contains 1 IQ domain.

## 画像



Detection limit for recombinant GST tagged IQCH is 1 ng/ml as a capture antibody.  
Antibody conc 10 µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors