

## Product datasheet

# Anti-Cdx1 antibody ab118358

画像数 3

### 製品の概要

製品名	Anti-Cdx1 antibody
製品の詳細	Mouse monoclonal to Cdx1
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB, ICC/IF, Sandwich ELISA
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment corresponding to amino acids 126-216 of Human Cdx1 with proprietary 26 kDa tag (NP_001795). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	HeLa cells

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Constituent: 99% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG2b
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab118358** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 28 kDa.

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 10 µg/ml.
Sandwich ELISA		Use a concentration of 10 µg/ml.

## ターゲット情報

機能	Could play a role in the terminal differentiation of the intestine.
組織特異性	Intestinal epithelium.
配列類似性	Belongs to the Caudal homeobox family. Contains 1 homeobox DNA-binding domain.
細胞内局在	Nucleus.


## 画像




Anti-Cdx1 antibody (ab118358) at 5 µg/ml +  
Immunizing peptide at 0.2 µg

**Predicted band size:** 28 kDa

The observed band size may not correspond to the predicted protein molecular weight as the immunogen (recombinant fragment) was used for the test lane.



ab118358, at 10 µg/ml staining Cdx1 in HeLa cells by Immunofluorescence.



Detection limit for recombinant tagged Cdx1 is 1 ng/ml as a capture antibody.

Sandwich ELISA - Anti-Cdx1 antibody (ab118358)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors