

## Product datasheet

# Anti-TBC1D21 antibody [7B1] ab119060

### 画像数 1

#### 製品の概要

製品名	Anti-TBC1D21 antibody [7B1]
製品の詳細	Mouse monoclonal [7B1] to TBC1D21
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Human
免疫原	Recombinant full length Human TBC1D21 produced in HEK293T cells (NP_699187).
ポジティブ・コントロール	HEK293T cell lysate transfected with pCMV6-ENTRY TBC1D21 cDNA
特記事項	Dilute in PBS (pH7.3) before use. Stable for 12 months from date of receipt.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 50% Glycerol, 1% BSA
精製度	Protein G purified
特記事項(精製)	ab119060 was purified from mouse ascites fluids by affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	7B1
アイソタイプ	IgG1

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab119060** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/2000. Predicted molecular weight: 39 kDa.

## ターゲット情報

機能	May act as a GTPase-activating protein for Rab family protein(s).
配列類似性	Contains 1 Rab-GAP TBC domain.

## 画像



**All lanes :** Anti-TBC1D21 antibody [7B1] (ab119060) at 1/2000 dilution

**Lane 1 :** HEK293T cell lysate transfected with pCMV6-ENTRY control

**Lane 2 :** HEK293T cell lysate transfected with pCMV6-ENTRY TBC1D21 cDNA

Lysates/proteins at 5 µg per lane.

**Predicted band size :** 39 kDa

HEK293T cell lysates were generated from transient transfection of the cDNA clone (RC202410)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors