

Product datasheet

Anti-C18orf1 antibody ab122994

画像数 1

製品の概要

製品名	Anti-C18orf1 antibody
製品の詳細	Rabbit polyclonal to C18orf1
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Chicken
免疫原	Synthetic peptide corresponding to a region within C terminal amino acids 181-230 (EVMGHHPGAS FLHHQRSNAH RGSRLQFQQN NAESTVPIK GKDRKPGNLV) of Human C18orf1 (NP_852148). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	U937 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituents: 97% PBS, 2% Sucrose
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122994** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
----------	-----------	------

WB	Use a concentration of 1 µg/ml. Predicted molecular weight: 26 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	--

ターゲット情報

配列類似性	Belongs to the PMEPA1 family. Contains 1 LDL-receptor class A domain.
細胞内局在	Cell membrane.

画像



Anti-C18orf1 antibody (ab122994) at 1 µg/ml
+ U937 cell lysate at 10 µg

Predicted band size : 26 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors