

Product datasheet

Anti-gtf3c3 antibody ab123053

画像数 1

製品の概要

| | |
|--------------|---|
| 製品名 | Anti-gtf3c3 antibody |
| 製品の詳細 | Rabbit polyclonal to gtf3c3 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Rat 交差が予測される動物種: Mouse, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Human, Pig |
| 免疫原 | Synthetic peptide corresponding to a region within C terminal amino acids 821-870 (LGLTHLAIHY YQKALALPPI VVEGIEVDQL DLRRDIA YNM) of Rat gtf3c3 (NP_001101709). Run BLAST with ExPASy Run BLAST with NCBI |
| ポジティブ・コントロール | Rat liver lysate. |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Constituents: 97% PBS, 2% Sucrose |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab123053** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 101 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

ターゲット情報

| | |
|-------|--|
| 機能 | Involved in RNA polymerase III-mediated transcription. Integral, tightly associated component of the DNA-binding TFIIIC2 subcomplex that directly binds tRNA and virus-associated RNA promoters. |
| 配列類似性 | Contains 11 TPR repeats. |
| 細胞内局在 | Nucleus. |

画像



Anti-gtf3c3 antibody (ab123053) at 1 µg/ml +
Rat liver lysate at 10 µg

Predicted band size: 101 kDa

Gel concentration: 6-18%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors