

Product datasheet

Anti-GRAP2 antibody ab123702

画像数 1

製品の概要

製品名	Anti-GRAP2 antibody
製品の詳細	Rabbit polyclonal to GRAP2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 130-159 of Human GRAP2 (NP_004801.1)
ポジティブ・コントロール	MCF-7 cell line lysates

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.
バッファー	Preservative: 0.09% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
特記事項(精製)	ab123702 was purified through a protein A column, followed by peptide affinity purification.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab123702** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/100 - 1/500. Predicted molecular weight: 38 kDa.

<b>機能</b>	Interacts with SLP-76 to regulate NF-AT activation. Binds to tyrosine-phosphorylated shc.
<b>配列類似性</b>	Belongs to the GRB2/sem-5/DRK family. Contains 1 SH2 domain. Contains 2 SH3 domains.

## 画像



Anti-GRAP2 antibody (ab123702) at 1/100 dilution + MCF-7 cell line lysates at 35 µg

**Predicted band size:** 38 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors