

Product datasheet

Anti-Osteopontin antibody ab14176

製品の概要

製品名	Anti-Osteopontin antibody
製品の詳細	Rabbit polyclonal to Osteopontin
由来種	Rabbit
アプリケーション	適用あり: ELISA
種交差性	交差種: Human
免疫原	Synthetic peptide: VYGLRSKSKKFRRPDI , corresponding to amino acids 164-179 of Human Osteopontin.
	Run BLAST with Run BLAST with

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole Serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab14176** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/5000.

ターゲット情報

機能	<p>Binds tightly to hydroxyapatite. Appears to form an integral part of the mineralized matrix. Probably important to cell-matrix interaction.</p> <p>Acts as a cytokine involved in enhancing production of interferon-gamma and interleukin-12 and reducing production of interleukin-10 and is essential in the pathway that leads to type I immunity.</p>
組織特異性	Bone. Found in plasma.
配列類似性	Belongs to the osteopontin family.
翻訳後修飾	<p>Extensively phosphorylated on clustered serine residues.</p> <p>N- and O-glycosylated.</p> <p>Phosphorylation sites are present in the extracellular medium.</p>
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors