

Product datasheet

Anti-CD5 antibody [4C7], prediluted ab17132

画像数 1

製品の概要

製品名	Anti-CD5 antibody [4C7], prediluted
製品の詳細	Mouse monoclonal [4C7] to CD5, prediluted
アプリケーション	適用あり: IHC-P, IHC-Fr
種交差性	交差種: Human
免疫原	Recombinant fusion protein corresponding to the external domain of Human CD5 molecule.
ポジティブ・コントロール	Tonsil

製品の特性

製品の状態	Prediluted
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.05% Sodium Azide Constituents: Inert stabilizer
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	4C7
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab17132** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		
IHC-Fr		

追加情報

IHC-P: Prediluted and ready to use. We suggest an incubation period of 30 minutes at room temperature. Perform heat mediated antigen retrieval with 10 mM citrate buffer, pH 6.0 before

commencing with IHC staining protocol.

IHC-Fr: Prediluted and ready to use. We suggest an incubation period of 30 minutes at room temperature.


Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能	May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.
配列類似性	Contains 3 SRCR domains.
細胞内局在	Membrane.

Anti-CD5 antibody [4C7], prediluted 画像



Immunohistochemical analysis of formalin-fixed paraffin-embedded human tonsil tissue sections, labelling CD5 with ab17132.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD5 antibody [4C7], prediluted (ab17132)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors