

## Product datasheet

# Anti-IL9 antibody ab18066

### 製品の概要

製品名	Anti-IL9 antibody
製品の詳細	Rabbit polyclonal to IL9
由来種	Rabbit
アプリケーション	<b>適用あり:</b> Neutralising, ELISA, WB
種交差性	<b>交差種:</b> Mouse
免疫原	Highly pure (>98%) recombinant full length protein (Mouse).

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: None Constituents: PBS, pH 7.4 (when reconstituted)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab18066** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Neutralising		
ELISA		
WB		

### 追加情報

ELISA: Use at a concentration of 0.5 µg/ml. antibody. In conjunction with compatible secondary reagents, allows the detection of 0.2 - 0.4 ng/well of recombinant murine

IL9.

Neut: Use at a concentration of 0.006 - 0.012 µg/ml to yield one-half maximal inhibition [ND50] of the biological activity of murine IL9 (0.30 ng/ml).

WB: Use at a concentration of 0.1 - 0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant murine IL9 is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions. Predicted molecular weight: 16 kDa.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

---

機能	Supports IL-2 independent and IL-4 independent growth of helper T-cells.
配列類似性	Belongs to the IL-7/IL-9 family.
細胞内局在	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors