

Product datasheet

Anti-Procaïnamide antibody [706] ab20768

製品の概要

製品名	Anti-Procaïnamide antibody [706]
製品の詳細	Mouse monoclonal [706] to Procaïnamide
由来種	Mouse
アプリケーション	適用あり: ELISA
免疫原	Chemical/ Small Molecule (Procaïnamide conjugated to KLH).
特記事項	The concentration of this product is lot specific

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: PBS, pH 7.4
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	706
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab20768** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

ELISA

追加情報

ELISA: May be used to identify and quantitate procaïnamide. Dilute in PBS or medium which is identical to that used in the assay system. Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性

Procainamide is used to correct irregular heartbeats to a normal rhythm and to slow an overactive heart. This allows the heart to work more efficiently. Procainamide produces its beneficial effects by slowing nerve impulses in the heart and reducing sensitivity of heart tissues.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors