

Product datasheet

Anti-Digoxin antibody [3301] ab20814

★★★★☆ 1 Abreviews [画像数 1](#)

製品の概要

製品名	Anti-Digoxin antibody [3301]
製品の詳細	Mouse monoclonal [3301] to Digoxin
由来種	Mouse
特異性	ab20814 shows 100% reactivity with digoxin and digitoxin
アプリケーション	適用あり: ELISA, RIA
種交差性	Not applicable
免疫原	Digoxin-BSA

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 150mM Sodium chloride, 15mM Potassium phosphate. pH 7.2
特記事項 (精製)	Purified from ascites by DEAE-chromatography.
ポリ/モノ	モノクローナル
クローン名	3301
ミエローマ	P3x63-Ag8.653
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab20814** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ELISA	★★★★☆	Use at an assay dependent dilution.
RIA		Use at an assay dependent dilution.

ターゲット情報

関連性

Digoxin is a beta blocker drug originally derived from the foxglove plant, *Digitalis lanata*. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).

画像



ELISA analysis of Human blood serum treated with 0 or 5 ng/ml Digoxin. Samples were incubated with ab20814 (1/90 in diluent), or control anti-Digoxin IgG for 2 hours at 20°C.

ELISA - Anti-Digoxin antibody [3301] (ab20814)

This image is courtesy of an Abreview submitted by Dr Julie Richards

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors