

Product datasheet

Anti-Verotoxin I/SLT 1b antibody [357] ab21093

製品の概要

製品名	Anti-Verotoxin I/SLT 1b antibody [357]
製品の詳細	Mouse monoclonal [357] to Verotoxin I/SLT 1b
由来種	Mouse
特異性	Specific for Verotoxin I - B subunit (SLT-1b).
アプリケーション	適用あり: Neutralising, ELISA
免疫原	Purified toxoid.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 10mM PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	357
ミエローマ	unknown
アイソタイプ	IgG1
軽鎖の種類	unknown

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab21093** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Neutralising		Use at an assay dependent dilution.
ELISA		1/20 - 1/200.

ターゲット情報

関連性	Shiga toxin contains a single A-chain and five copies of a B-chain. The B-chain is responsible for binding the holotoxin to specific receptors on the target cell surface.
細胞内局在	Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors