

## Product datasheet

### Anti-AU5 tag antibody [AU5] ab24576

#### 1 Abreviews

#### 製品の概要

製品名	Anti-AU5 tag antibody [AU5]
製品の詳細	Mouse monoclonal [AU5] to AU5 tag
由来種	Mouse
アプリケーション	<b>適用あり:</b> WB, ICC/IF, IP
免疫原	This antibody was raised against an SDS disrupted BPV1 preparation and mapped to the six amino acid epitope TDFYLK.
エピトープ	This antibody maps to the six amino acid epitope TDFYLK.
特記事項	<p>This clone was specifically designed for use in epitope tagging. Because of its small size, clone AU5 is unlikely to alter the activity of the cloned sequence.</p> <p>Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. If you would like help finding an alternative product, please do not hesitate to contact our scientific support team.</p>

#### 法規制情報

**医薬用外毒物**

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.03% Thimerosal (merthiolate) Constituent: PBS
精製度	Protein G purified
一次抗体 備考	This clone was specifically designed for use in epitope tagging. Because of its small size, clone AU5 is unlikely to alter the activity of the cloned sequence.
ポリ/モノ	モノクローナル
クローン名	AU5

## アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab24576** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500. Predicted molecular weight: 1 kDa.
ICC/IF		1/500.
IP		1/150.

## ターゲット情報

**関連性** AU5 is a commonly used epitope tag engineered onto the N- or C- terminus of a protein of interest so that the tagged protein can be analyzed and visualized. Because of the small size of the epitope, it is unlikely to alter the activity of the cloned sequence.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors