

## Product datasheet

# Anti-NKG2A + C+ E antibody [20d5] (FITC) ab24866

### 製品の概要

製品名	Anti-NKG2A + C+ E antibody [20d5] (FITC)
製品の詳細	Rat monoclonal [20d5] to NKG2A + C+ E (FITC)
由来種	Rat
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	Tissue/ cell preparation: CHO cells transfected with murine NKG2.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium azide Constituent: PBS
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	20d5
アイソタイプ	IgG2a
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab24866** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

### 追加情報

Flow Cyt: Use 1µg for 10<sup>6</sup> cells.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## ターゲット情報

---

関連性	NKG2A, NKG2C and NKG2E, are isoforms of the NKG2 lectin-like family. In mice, NKG2 subunits associate with CD94 to form heterodimers at the surface of natural killer (NK) cells. The CD94/NKG2 heterodimer is the receptor for a non-classical MHC class I ligand, which is Qa 1 in the mouse. NKG2 is expressed on NK cells, lymphokine-activated killer (LAK) T cells, and some CD8+ memory T cells.
細胞内局在	Type II membrane protein.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors