

Product datasheet

Anti-Ia Antigen antibody [Cla] (Biotin) ab25158

1 References

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Ia Antigen antibody [Cla] (Biotin) |
| 製品の詳細 | Mouse monoclonal [Cla] to Ia Antigen (Biotin) |
| 標識 | Biotin |
| 特異性 | ab25158 recognises Ia antigen (65 kDa glycoprotein consisting of a 32 kDa α -chain non-covalently bound to a 27 kDa β -chain). |
| アプリケーション | 適用あり: Flow Cyt, IP, Inhibition Assay |
| 種交差性 | 交差種: Chicken |
| 免疫原 | The details of the immunogen for this antibody are not available. |
| 特記事項 | The molecule recognized by this antibody is similar in structure to the murine Ia antigens and human DR antigens. |

製品の特性

| | |
|---------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Preservative: 0.1% Sodium Azide Constituents: PBS |
| 精製度 | Purified IgM |
| 一次抗体 備考 | The molecule recognized by this antibody is similar in structure to the murine Ia antigens and human DR antigens. |
| ポリ/モノ | モノクローナル |
| クローン名 | Cla |
| アイソタイプ | IgM |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25158** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|------------------|-----------|---|
| Flow Cyt | | Use 1µg for 10 ⁶ cells. ab18396 -Mouse monoclonal IgM, is suitable for use as an isotype control with this antibody. |
| IP | | Use at an assay dependent concentration. |
| Inhibition Assay | | Use at an assay dependent concentration. Inhibition of MLR. |

ターゲット情報

| | |
|-------|---|
| 関連性 | The chicken Ia antigen is expressed on B cells, a subpopulation of the monocyte/macrophage lineage of cells, and mitogen-activated T cells. |
| 細胞内局在 | Type I membrane protein. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors