

## Product datasheet

# Anti-Qa 2 antibody [695H1-9-9] (Biotin) ab25236

## 2 References

### 製品の概要

製品名	Anti-Qa 2 antibody [695H1-9-9] (Biotin)
製品の詳細	Mouse monoclonal [695H1-9-9] to Qa 2 (Biotin)
由来種	Mouse
標識	Biotin
アプリケーション	<b>適用あり:</b> Flow Cyt, Functional Studies
種交差性	<b>交差種:</b> Mouse
免疫原	Tissue/ cell preparation: C3H.SW mouse skin graft and splenocytes.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	695H1-9-9
アイソタイプ	IgG2a
軽鎖の種類	kappa

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25236** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		

Functional Studies

#### 追加情報

Flow Cyt: Use 1µg for 10<sup>6</sup> cells.

Complement mediated cytotoxicity experiments: Use at an assay dependent dilution.

Modulation of Qa2 expression from peripheral T lymphocytes: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

#### ターゲット情報

#### 関連性

Qa antigens are 'non classical' MHC class I molecules of mice. Qa2 is present on restricted functional lymphocyte subpopulations. The protein product of the preimplantation embryo development (Ped) gene, Qa2 influences the rate of embryonic cleavage during the preimplantation stages of development.

#### 細胞内局在

Plasma membrane

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors