

Product datasheet

Anti-CD5 antibody [f43] ab25337

製品の概要

| | |
|----------|---|
| 製品名 | Anti-CD5 antibody [f43] |
| 製品の詳細 | Mouse monoclonal [f43] to CD5 |
| 由来種 | Mouse |
| 特異性 | ab25337 recognises feline CD5. |
| アプリケーション | 適用あり: Flow Cyt, IHC-Fr, IP |
| 種交差性 | 交差種: Cat |
| 免疫原 | The details of the immunogen for this antibody are not available. |
| 特記事項 | This antibody has been shown to be useful as a pan T cell marker. |

製品の特性

| | |
|---------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | pH: 8.20 Constituent: 100% Borate buffered saline |
| 精製度 | IgG fraction |
| 一次抗体 備考 | This antibody has been shown to be useful as a pan T cell marker. |
| ポリ/モノ | モノクローナル |
| クローン名 | f43 |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25337** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| Flow Cyt | | Use at an assay dependent concentration. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody. |
| IHC-Fr | | Use at an assay dependent concentration. Acetone fixed. |
| IP | | Use at an assay dependent concentration. |

ターゲット情報

| | |
|-------|--|
| 機能 | May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2. |
| 配列類似性 | Contains 3 SRCR domains. |
| 細胞内局在 | Membrane. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors