

Product datasheet

Anti-CD5 antibody [f43] ab25337

製品の概要

製品名	Anti-CD5 antibody [f43]
製品の詳細	Mouse monoclonal [f43] to CD5
特異性	ab25337 recognises feline CD5.
アプリケーション	適用あり: Flow Cyt, IHC-Fr, IP
種交差性	交差種: Cat
免疫原	The details of the immunogen for this antibody are not available.
特記事項	This antibody has been shown to be useful as a pan T cell marker.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 8.20 Constituent: 100% Borate buffered saline
精製度	IgG fraction
一次抗体 備考	This antibody has been shown to be useful as a pan T cell marker.
ポリ/モノ	モノクローナル
クローン名	f43
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25337** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. ab170190 -Mouse monoclonal IgG1, is suitable for use as an

isotype control with this antibody.

IHC-Fr	Use at an assay dependent concentration. Acetone fixed.
--------	------------------------------------------------------------

IP	Use at an assay dependent concentration.
----	------------------------------------------

ターゲット情報

機能	May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.
配列類似性	Contains 3 SRCR domains.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors