

Product datasheet

Anti-LAMP1 antibody [H4A3] (Phycoerythrin) ab25367

1 References [画像数 1](#)

製品の概要

製品名	Anti-LAMP1 antibody [H4A3] (Phycoerythrin)
製品の詳細	Mouse monoclonal [H4A3] to LAMP1 (Phycoerythrin)
由来種	Mouse
標識	Phycoerythrin. Ex: 488nm, Em: 575nm
アプリケーション	適用あり: Flow Cyt, IHC-Fr
種交差性	交差種: Human
免疫原	The details of the immunogen for this antibody are not available.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 16% Sucrose, PBS; stabilizing agent
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	H4A3
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25367** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt	Use 10µl for 10 ⁶ cells. ab91357 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
----------	---


IHC-Fr

Use at an assay dependent concentration.

ターゲット情報

機能	Presents carbohydrate ligands to selectins. Also implicated in tumor cell metastasis.
配列類似性	Belongs to the LAMP family.
翻訳後修飾	O- and N-glycosylated; some of the 18 N-linked glycans are polylactosaminoglycans.
細胞内局在	Cell membrane. Endosome membrane. Lysosome membrane. This protein shuttles between lysosomes, endosomes, and the plasma membrane.

画像



Flow cytometry analysis staining LAMP1 in Jurkat cells using ab25367 at a dilution of 10 uL/10⁶ cells. Mouse anti-human CD107a-PE was used as an isotype control.

Flow Cytometry - Anti-LAMP1 antibody [H4A3]
(Phycoerythrin) (ab25367)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors