

Product datasheet

Anti-CD90 / Thy1 antibody [30-H12] (Phycoerythrin) ab25694

製品の概要

| | |
|----------|--|
| 製品名 | Anti-CD90 / Thy1 antibody [30-H12] (Phycoerythrin) |
| 製品の詳細 | Rat monoclonal [30-H12] to CD90 / Thy1 (Phycoerythrin) |
| 由来種 | Rat |
| 標識 | Phycoerythrin. Ex: 488nm, Em: 575nm |
| アプリケーション | 適用あり: Flow Cyt, IHC-Fr |
| 種交差性 | 交差種: Mouse |
| 免疫原 | Tissue, cells or virus corresponding to Mouse CD90/ Thy1. Immunogen is Mouse thymus or spleen. |

特記事項

This antibody has been shown to be useful in 'in vivo' and 'in vitro' studies of T cell depletion.

製品の特性

| | |
|---------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | pH: 7.2 Preservative: 0.09% Sodium azide Constituents: 0.16% Sodium phosphate, 0.87% Sodium chloride, 0.1% Gelatin |
| 精製度 | Affinity purified |
| 一次抗体 備考 | This antibody has been shown to be useful in 'in vivo' and 'in vitro' studies of T cell depletion. |
| ポリ/モノ | モノクローナル |
| クローン名 | 30-H12 |
| アイソタイプ | IgG2b |
| 軽鎖の種類 | kappa |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25694** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| Flow Cyt | | Use 0.05µg for 10 ⁶ cells. |
| IHC-Fr | | Use at an assay dependent concentration. |

ターゲット情報

| | |
|-------|---|
| 機能 | May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain. |
| 配列類似性 | Contains 1 Ig-like V-type (immunoglobulin-like) domain. |
| 細胞内局在 | Cell membrane. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors