

Product datasheet

Anti-Caspase 5 antibody ab25899

画像数 1

製品の概要

製品名	Anti-Caspase 5 antibody
製品の詳細	Rabbit polyclonal to Caspase 5
由来種	Rabbit
特異性	ab25899 recognises Caspase-5.
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide, consisting of 16 N-terminal amino acids of Human Caspase 5
ポジティブ・コントロール	Ramos cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.02% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab25899** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

WB

追加情報

WB: Use at a concentration of 1 µg/ml. Detects a band of approximately 48 kDa (predicted molecular weight: 48 kDa).

Depending on cell lines or tissues used, other cleavage products may be observed.

Can be blocked with [Caspase 5 peptide \(ab39726\)](#).


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

機能	Mediator of programmed cell death (apoptosis).
組織特異性	Expressed in barely detectable amounts in most tissues except brain, highest levels being found in lung, liver and skeletal muscle.
配列類似性	Belongs to the peptidase C14A family. Contains 1 CARD domain.
翻訳後修飾	The two subunits are derived from the precursor sequence by an autocatalytic mechanism.

画像



Lane 1 : Anti-Caspase 5 antibody (ab25899)
at 0.5 µg/ml

Lane 2 : Anti-Caspase 5 antibody (ab25899)
at 1 µg/ml

Lane 3 : Anti-Caspase 5 antibody (ab25899)
at 2 µg/ml

All lanes : Ramos cells

Predicted band size: 48 kDa

Observed band size: 48 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors