

Product datasheet

Anti-Transglutaminase 5 antibody ab26992

1 References

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Transglutaminase 5 antibody |
| 製品の詳細 | Rabbit polyclonal to Transglutaminase 5 |
| 由来種 | Rabbit |
| 特異性 | This antibody is specific to Transglutaminase 5 protein from keratinocyte cells. |
| アプリケーション | 適用あり: WB |
| 種交差性 | 交差種: Mouse, Human |
| 免疫原 | A 15 residue synthetic peptide derived from a domain specific for Transglutaminase 5 protein. |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| バッファー | Preservative: None Constituents: Whole serum |
| 精製度 | Whole antiserum |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab26992** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| WB | | 1/200 - 1/1000. Predicted molecular weight: 80 kDa. |

ターゲット情報

| | |
|--------------|--|
| 機能 | Catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins. Contributes to the formation of the cornified cell envelope of keratinocytes. |
| 組織特異性 | Expressed in foreskin keratinocytes. |
| 関連疾患 | Defects in TGM5 are a cause of peeling skin syndrome type A (APSS) [MIM:609796]. A non-inflammatory form of peeling skin syndrome, a genodermatosis characterized by the continuous shedding of the outer layers of the epidermis. In APSS patients, skin peeling is strictly limited to the dorsa of the hands and feet, and it is accompanied by painless erythema and spontaneous non-scarring healing. Ultrastructural and histological analysis shows a level of blistering high in the epidermis at the stratum granulosum-stratum corneum junction. |
| 配列類似性 | Belongs to the transglutaminase superfamily. Transglutaminase family. |
| 細胞内局在 | Cytoplasm. Associated with intermediate filaments. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors