

Product datasheet

Anti-FLJ12644 antibody ab28728

画像数 1

製品の概要

製品名	Anti-FLJ12644 antibody
製品の詳細	Rabbit polyclonal to FLJ12644
由来種	Rabbit
特異性	ab28728 recognises FLJ12644.
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human
免疫原	A region within synthetic peptide: ECEKAYFYMS CLVKHKRIHS REKRGDSVKV ENPSTASHSL SPSEHVQGKS, corresponding to C terminal amino acids 434-483 of Human FLJ12644 Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Jurkat Cell Lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab28728** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ELISA		1/312500.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1.25 µg/ml. Predicted molecular weight: 58 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Transcriptional repressor. Regulator of transcriptional factor complexes and may suppress SRE and AP-1 transcription activities mediated by growth factor signaling pathways.
組織特異性	Highly expressed in heart, skeletal muscle, and brain. Lower expression in liver, lung, kidney, pancreas and placenta.
配列類似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 10 C2H2-type zinc fingers. Contains 1 KRAB domain.
発生段階	Expressed in fetal heart, brain, placenta, lung, liver, skeletal muscle, kidney, and pancreas.
ドメイン	The KRAB domain is required for transcriptional repression.
細胞内局在	Nucleus.

画像



Anti-FLJ12644 antibody (ab28728) at 1.25 µg/ml + Jurkat cell lysate

Secondary

HRP conjugated anti-rabbit IgG at 1/50000 dilution

Predicted band size: 58 kDa

Observed band size: 60 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors