

## Product datasheet

# Anti-DDX27 antibody ab36312

画像数 1

### 製品の概要

製品名	Anti-DDX27 antibody
製品の詳細	Rabbit polyclonal to DDX27
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Cow, Dog
免疫原	A region within synthetic peptide: KKKRAATTLD EKIEKVRKKR KTEDKEAKSG KLEKEKEAKE GSEPKQEDL, corresponding to N terminal amino acids 109-158 of Human DDX27 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	HepG2 cell lysate.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Protein A purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab36312** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/1.5625e+006.


アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

WB Use a concentration of 0.625 µg/ml. Detects a band of approximately 60 kDa (predicted molecular weight: 88 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## ターゲット情報

機能	Probable ATP-dependent RNA helicase.
配列類似性	Belongs to the DEAD box helicase family, DDX27/DRS1 subfamily. Contains 1 helicase ATP-binding domain. Contains 1 helicase C-terminal domain.
細胞内局在	Nucleus.

## 画像



Anti-DDX27 antibody (ab36312) at 0.625 µg/ml + HepG2 cell lysate at 10 µg

### Secondary

HRP conjugated anti Rabbit IgG 1/50,000 - 1/100,000

**Predicted band size** : 88 kDa

**Observed band size** : 60 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors