


## Product datasheet

# Anti-PARP11 antibody ab42485

### 画像数 1

#### 製品の概要

製品名	Anti-PARP11 antibody
製品の詳細	Rabbit polyclonal to PARP11
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Cow 
免疫原	Synthetic peptide (Human). More details of this immunogen are not available.
ポジティブ・コントロール	Fetal muscle lysate. Fetal lung lysate.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	Peptide affinity chromatography.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab42485** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 39 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


## ターゲット情報

---

**関連性** Poly(ADP-ribosyl)ation is a DNA-damage-dependent post-translational modification of histones. Poly(ADP-ribose) polymerases (PARPs) are a family of 18 proteins, encoded by different genes and displaying a conserved catalytic domain.

## 画像

---



**All lanes :** Anti-PARP11 antibody (ab42485)  
at 1 µg/ml

**Lane 1 :** Fetal muscle lysate

**Lane 2 :** Fetal lung lysate

**Predicted band size :** 39 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors