

Product datasheet

Anti-Znf206 antibody ab45344

2 Abreviews

製品の概要

製品名 Anti-Znf206 antibody
 製品の詳細 Rabbit polyclonal to Znf206

i This product is a **fast track antibody**. It has been affinity purified and shows high titre values against the immunizing peptide by ELISA.

[Read the terms of use »](#)

由来種 Rabbit

種交差性

交差が予測される動物種: Mouse, Human **!**

免疫原 Synthetic peptide conjugated to KLH derived from within residues 600 - 700 of Human Znf206.Immunogen の所有権に関して(Peptide available as [ab45343](#).)

製品の特性

製品の状態 Liquid
 保存方法 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
 バッファー Preservative: 0.02% Sodium Azide
 Constituents: 1% BSA, PBS, pH 7.4
 精製度 Immunogen affinity purified
 ポリ/モノ ポリクローナル
 アイソタイプ IgG

アプリケーション

追加情報 This antibody gave a positive result in ELISA against the immunizing peptide ([ab45343](#)).
 Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性	Znf206 is a zinc finger protein and probable transcriptional regulator. It is expressed in embryonic stem cells and downregulated upon differentiation. Both Oct4 and Sox2 bind the znf206 gene.
細胞内局在	Nuclear

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors