

Product datasheet

Placenta (mouse) normal tissue slides ab4612

製品の概要

製品名	Placenta (mouse) normal tissue slides
特記事項	After excision, tissue was fixed in paraformaldehyde and processed to paraffin wax. Tissue sections were cut to a thickness of 5µm and mounted on positively charged glass slides. One tissue from each lot was stained by haematoxylin and eosin (H+E) and then microscopically examined to ensure quality. All tissue is examined / diagnosed by a licensed pathologist and is ethically obtained. This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.
アプリケーション	適用あり: IHC-P, In situ hybridization

製品の特性

保存方法	Store at +4°C.
スライド 備考	This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab4612** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent concentration.
In situ hybridization		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors