


Product datasheet

Anti-Endothelin 2 antibody ab47923

製品の概要

製品名	Anti-Endothelin 2 antibody
製品の詳細	Rabbit polyclonal to Endothelin 2
由来種	Rabbit
特異性	This antibody can cross react with Endothelin-1 (20%) and Endothelin-3 (less than 10%).
アプリケーション	適用あり: RIA, IP, ELISA
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Cow, Cat, Dog, Ferret 
免疫原	Synthetic peptide: CSCSSWLDKECVYFCHLDIIW conjugated to KLH via Cys, corresponding to amino acids 49-69 of Human Endothelin 2. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>

法規制情報

**医薬用外毒物**

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 50% Glycerol, PBS, pH 7.5
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab47923** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
RIA		Use at an assay dependent dilution.
IP		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.

#### ターゲット情報

機能	Endothelins are endothelium-derived vasoconstrictor peptides.
組織特異性	Expressed in lung, but not in placental stem villi vessels or cultured placental villi smooth muscle cells.
配列類似性	Belongs to the endothelin/sarafotoxin family.
細胞内局在	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors